

nPOD Online Pathology 101

Martha Campbell-Thompson

April 24, 2013

Welcome

- Asking questions:
 - Welcome at any step- others likely have the same
 - How- webinar chat function send to all-
 - If I don't read your question, talk
- Please mute your phones to avoid background noise- use mute button or *69


Best Viewing

- Browser: Microsoft Internet Explorer
 - Known issues with Chrome and Firefox
- For access: see [JDRF nPOD/For Investigators/Online Pathology Information](#)
 - Users limited to those in research
 - Internet security requires individual accounts with no sharing
- Webinar available on the webpage

Digital Pathology Primer

- Background: Case Processing
 - Histology and IHC
 - Pathology Review
- Spectrum: Information Management System
 - Cases- Specimens - Slides
 - Advanced searches
 - Export data
 - Take photos
 - Link H&E and IHC images

Case Processing



Pancreas: head, body, and tail

Paraffin blocks- 5-10 (OCT takes precedence over fixed when sample size is small)

OCT blocks- 5-10 all regions or as indicated by size

PLN:

OCT blocks – 3-5- as many as feasible depending on total numbers




Cells- place several dissected LNs in 15ml tubes containing sterile RPMI, hold at refrigeration until shipped or processed. 1-4 tubes depending on numbers

Frozen vials as for pancreas. Paraffin is optional.

Spleen and NonPLN: as for PLN.

Duodenal mucosa: as for pancreas.

Pancreas Subsections




Standard Histology and IHC

3-4 unstained serial Sections/block, 4um	Blocks 02A, 04A PanHead 02A, 04A PanBody, Tail 01 Spleen 01 Bone Marrow, Lesion(s)
H&E	1. H&E
IHC	2. Ki67+Insulin 3.CD3+Glucagon 4. PP (PanHead only)
Slide scanning	All stained slides

nPOD Slide Labels

Original Labels



Case ID
Block Number
Sample ID
Stain

Current Labels



Case ID = 6189
Block Number= 01
Sample ID= Spleen
Stain = H&E

Within 2D code:

Slide Number
Aliquot ID (paraffin or OCT)
Date

Case Information

CaseID:	6070	
Group:	T1D	"Pending" until lab data returned
AutoAb (Denver):	IA-2A+ mIAA+	Denver Autoantibody Core (G. Eisenbarth)
Age (years):	22.6 (7 w/diabetes)	
Gender:	Female	
Ethnicity:	Caucasian	
C-peptide (ng/ml):	<0.05	Northwest Lipid Research (S. Marcovina)
HbA1c:		
UNOS_BMI:	21.6	
ClinicalHistory:		"Accepted as" information, other
Histopathology:	Ins+ islets (few); occ. Ki67+ islet cell. CD3+ Insulitis. Alternating blocks can be devoid of beta cells.	Synopsis of findings starting with islets
HiRes HLA:	A*0201, 0205 B*3801, 5801 C*0701, 1203 DRB1*1001, 1601 DQA1*0101, 0102 DQB1*0501, 0502 DPA1*0103, 0201 DPB1*1301, 1401	Oakland Children's (J. Noble)
Transplant HLA:	A*02/02, B*38/58, DRB1*10/16, DQB1*05/05	From chart
UNOS HLA Class I :	A1*02; B1*38; Cw1*07; A2*No second antigen detected; B2*58; Cw2*07; Bw4+; Bw6-	
UNOS HLA Class II :	DR1*10; DQ1*05; DPw1*Not Tested; DR2*16; DQ2*No second antigen detected; DPw2*Not Tested; DR51+; DR52- ; DR53-;	
		UNOS- Data uploaded via Data Management Core (J. Kaddis)

Case Searches

Spectrum™ [Advanced Search](#)

Cases Specimens Digital Slides Administrative Help Log off

Load Saved Search (Optional)


Modify Search Criteria

Search For

The Case Fields

CaselD
Group
AutoAb (Denver)
Age (years)
Gender
Ethnicity
C-peptide (ng/ml)
HbA1c
UNOS_BMI
ClinicalHistory
Histopathology
HiRes HLA
Transplant HLA
UNOS HLA Class I

Save



Note - search names containing an asterisk (*) will prompt for modification before execution

Histopathology Searches

Load Saved Search (Optional)

-- Select --

Modify Search Criteria

Search For

+ The Case Fields

☐ Histopathology

+ The Specimen Fields

+ The Digital Slide Fields

Save Search (Optional)

Save Search Criteria as

Note - search names containing an asterisk (*) will prompt for modification before execution

Specimens = Tissues = SampleID

Case Specimens

[View Images](#) | [Open Data](#)

<input type="checkbox"/>	Icon	Image	ID ↓	SampleID	Subproject	Comment
<input type="checkbox"/>		n/a	1998	PanHead		

[View Images](#) | [Open Data](#)

<input type="checkbox"/>	Icon	Label	Image	Slide ID	Block ID	CaseID	Block #	Slide Number	SampleID	Stain	Aliquot ID	Date	Comment	Analys
<input type="checkbox"/>				51316		6189	02A	01	PanHead	H&E	Paraffin	3/20/12		
<input type="checkbox"/>				51317		6189	04A	01	PanHead	H&E	Paraffin	3/20/12		
<input type="checkbox"/>				51449		6189	02A	06	PanHead	No primary	Paraffin	3/22/12		
<input type="checkbox"/>				51364		6189	02A	02	PanHead	Ki67+Insulin	Paraffin	3/22/12		
<input type="checkbox"/>				51375		6189	04A	02	PanHead	Ki67+Insulin	Paraffin	3/22/12		
<input type="checkbox"/>				51354		6189	02A	05	PanHead	Pan Polypeptide	Paraffin	3/22/12		
<input type="checkbox"/>				51356		6189	02A	03	PanHead	CD3+Glucagon	Paraffin	3/22/12		
<input type="checkbox"/>				51357		6189	04A	03	PanHead	CD3+Glucagon	Paraffin	3/22/12		

50

<input type="checkbox"/>		n/a	1997	PanBody		
<input type="checkbox"/>		n/a	1996	PanTail		
<input type="checkbox"/>		n/a	1995	Spleen		

Slide Information

Modify Search Criteria

Search For

☐ The Case Fields

☐ The Specimen Fields

☐ The Digital Slide Fields

Slide ID
Block ID
CaseID
Block #
Slide_Number
SampleID
Stain
Aliquot ID
Date
Comment
Analysis Progress

Save

Save

Delete

an asterisk (*) will prompt for modification before execution

Specimen & Slide Searches

Load Saved Search (Optional)

-- Select --

Modify Search Criteria

Search For

+ The Case Fields

+ The Specimen Fields

☐ SampleID equal to

+ The Digital Slide Fields


☐ Comment contains

Save Search (Optional)


Save Search Criteria as

Note - search names containing an asterisk (*) will prompt for modification before execution

WebScope- Macs or Default PC



ImageScope



iPhone, iPad

Aperio® ePathViewer Enables Easy Access to eSlides

The Aperio ePathViewer application allows mobile access to eSlide images from stored pathology data sets and is designed for optimal use with the iPad® mobile digital device, iPad2® mobile digital device, and the new iPad® mobile digital device. A separate app is available for use on the iPhone® 4 mobile digital device and iPhone® 4S mobile digital device.

The application was carefully developed to provide important capabilities to support remote viewing for ePathology slides. Users can select images from active thumbnails of eSlide images, adjust the fixed magnification up to the maximum from the original scan, and display eSlide metadata. In addition, the current view snapshot can be emailed if needed.

For quick reference, the application has built-in default websites, including the [Juan Rosai Collection of Surgical Pathology Seminars](#), as well as Aperio sample images from College of American Pathologists (CAP).



Aperio ePathViewer home page displays a website panel. Moving through select websites and images is as easy as a finger swipe.

For ease of use, a single tap toggles the display of an image set, toolbar, image label, thumbnail and color balance. More instructions on use are available below.

Download the ePathViewer Instruction Guide Here.



First Name*:

Last Name*:

Organization*:

Job Title*:

Phone*:

Email*:

Postal Code*:

Which Best Describes Your Area of Interest?*

☐ Healthcare ☐ Life Science

Call Me*: ☐ Yes ☐ No

Where Did You Hear About Aperio?*

Read the ePathViewer FAQs

- Get the ePathViewer App from the App Store:
 - Download ePathViewer App for the iPad
 - Download ePathViewer App for the iPhone

Request Account nPOD Website

nPOD Newsletter

A quarterly e-newsletter for nPOD investigators, partners, and friends of the project. [View the nPOD Newsletter archive.](#)

► Sign Up Here

Contact

For organ placement and donor suitability, contact:

1-866-731-6585

We are available 24/7.

For other inquiries, please contact:
(352) 273-8277
nPOD@pathology.ufl.edu

► Additional Contact Info

nPOD Online Pathology Investigator Request Form

Name:	*	<input type="text"/>
Phone Number:	*	<input type="text"/>
Email Address:	*	<input type="text"/>
Confirm Email Address:	*	<input type="text"/>
Name of Institution:	*	<input type="text"/>
Reason for Request:	*	<input type="text"/>
		<input type="button" value="Submit"/> <input type="button" value="Cancel"/>

Mingder Yang return your user information via email.

Customize Your Account

The screenshot shows a web browser window with the address bar displaying `http://ahc-path-aperio11.ad.ufl.edu/E`. The browser's address bar and tabs show the page title "Spectrum - User Settings". The browser's menu bar includes "File", "Edit", "View", "Favorites", "Tools", and "Help". The browser's toolbar shows various icons, including a search icon, a document icon, and a folder icon. The browser's status bar shows the following text: "par-12-017 Shared Instru...", "IDS", "ITN", "Campbell-Thompson, Ma...", "Spectrum", "UF MPC", "Reference", "OPPC", "nPOD", "Pathology", and "PATH". The main content area of the browser shows the "Spectrum" logo, a search bar with the text "-- Search --", and a link to "Advanced Search". The user's name, "User: Martha Campbell-Thompson", is displayed in the top right corner. Below the search bar, there is a navigation menu with the following items: "Cases", "Specimens", "Digital Slides", "Administrative", "Help", and "Log off". The "Administrative" item is currently selected. Below the navigation menu, there is a green banner that reads "User data successfully updated". Below the banner, there is a blue button labeled "My Settings".

User Details

Login:	thompmc
Authentication Type:	Spectrum
Full Name:	Martha Campbell-Thompson
Phone Number:	<input type="text"/>
E-mail Address:	<input type="text" value="mct@ufl.edu"/>
New Password:	<input type="password"/>
Retype New Password:	<input type="password"/>
Spectrum Start Page	<input type="text" value="Welcome Page (default)"/>
Enable Auto-View Images:	<input checked="" type="checkbox"/>
Clinical Viewing Mode:	<input type="checkbox"/>

Let's Go!

